



Email Website info@adipogen.com www.adipogen.com

PRODUCT DATA SHEET

PBS Endotoxin-free (sterile)

Cat. No.: IAX-900-001 **Date:** 08-Jan-2013

FORMULATION: 1.5ml (LD15) or 10ml (L010) PBS pH 7.2 in sterile filtered double-distilled water (ddWater), without any additives.

ENDOTOXIN-FREE: <0.0002EU/μg, after repeated purification cycles over activated charcoal.

APPLICATION: Use for the dilution of TLR*pure*TMLPS and solubilisation of endotoxin-free reagents (e.g. ODNs).

HANDLING: Keep sterile.

SHIPPING: Ambient. STORAGE: 4°C.

STABILITY: 2 years after receipt.