

## PRODUCT DATA SHEET

### ddWater Endotoxin-free (sterile)

Cat. No.: IAX-900-002

Date: 08-Jan-2013

<b>FORMULATION:</b>	1.5ml (LD15) or 10ml (L010) sterile filtered double-distilled water (ddWater), without any additives.
<b>ENDOTOXIN-FREE:</b>	<0.0002EU/μg, after repeated purification cycles over activated charcoal.
<b>APPLICATION:</b>	Use for the dilution of TLR <i>pure</i> <sup>TM</sup> LPS, MPLA and Lipid A and solubilisation of endotoxin-free reagents (e.g. ODNs).
<b>HANDLING:</b>	Keep sterile.
<b>SHIPPING:</b>	Ambient.
<b>STORAGE:</b>	4°C.
<b>STABILITY:</b>	2 years after receipt.

**DISCLAIMER:** THIS PRODUCT IS NOT INTENDED OR APPROVED FOR HUMAN, DIAGNOSTICS OR VETERINARY USE. USE OF THIS PRODUCT FOR HUMAN OR ANIMAL TESTING MAY BE EXTREMELY HAZARDOUS AND MAY RESULT IN DISEASE, SEVERE INJURY, OR DEATH. THIS PRODUCT IS FOR RESEARCH USE ONLY (RUO).

**MATERIAL SAFETY DATA:** This material should be considered hazardous until information to the contrary becomes available. Do not ingest, swallow, inhale or get into the blood stream. Do not get in eyes, on skin, or clothing. Wash thoroughly after handling. This information contains some, but not all, of the information required for the safe and proper use of this material. Access to this material must be restricted to personnel, who is appropriately experienced, qualified, competent and properly trained to use it. Material Safety Data Sheet is available upon request.