



## **PRODUCT DATA SHEET**

Physiological Saline [Sodium Chloride 0.9%] Endotoxin-free (sterile)

**Cat. No.:** IAX-900-003 **Date:** 10-October-2016

FORMULATION: 1.5ml (LD15) or 10ml (L010) NaCl pH 5.5-6.0 in sterile double-distilled water (ddWater), without any additives.

**ENDOTOXIN-FREE:** <0.0002EU/µg, after repeated purification cycles over activated charcoal.

**APPLICATION:** Use for the dilution of TLR*pure*<sup>TM</sup>LPS and solubilisation of endotoxin-free reagents (e.g. ODNs).

**HANDLING:** Keep sterile.

SHIPPING: Ambient. STORAGE: 4°C.

**STABILITY:** 2 years after receipt.